

Sniffing Out the Standards: Making Sense of Evolving Nasal Spray Testing

From Regulatory Guidance to Real-World Practice

Maria Smith, Direction of Applications & Business Development
Proveris Laboratories



IPAC-RS Event
Friday, September 19th

Agenda

- The Growing Importance of Nasal Delivery
- Regulatory Landscape: Global & Evolving Standards
- Making Sense of the Standards
(What They Really Mean in Practice)
- Case Studies & Special Considerations
(PBE vs. ABE, Device Testing, Actuation)
- Ensuring Success: From Method Validation to Spec Setting
- Looking Forward: Pain Points, Unmet Needs &
Next-Gen Testing Approaches
- Key Takeaways

The Growing Importance of Nasal Delivery



Growing importance of nasal drug delivery in therapeutic innovation - anticipated \$12B+ growth over the next 10 years¹



Increased prevalence of respiratory disease and shift toward non-invasive drug delivery methods



North America and EU hold majority market share



Regulatory standards and expectations continue to evolve as technological advancements are made, and innovation solutions are developed

¹<https://www.futuremarketinsights.com/reports/nasal-spray-market>

Regulatory Overview

U.S FDA

- *Spray and Inhalation Solution, Suspension, and Spray Drug Products-Chemistry, Manufacturing, and Controls Documentation - July 2002*
- *Bioavailability and Bioequivalence Studies for Nasal Aerosols and Nasal Sprays for Local Action -April 2003 DRAFT*
- *48 Product Specific Guidance - All in DRAFT*

USP General Chapters

- *USP <5> Inhalation and Nasal Drug Products—General Information and Product Quality Tests - Aug 2023*
- *USP <601> Inhalation and Nasal Drug Products: Aerosols, Sprays, and Powders—Performance Quality Tests - May 2024*

EMA

- *Guideline on the pharmaceutical quality of inhalation and nasal medicinal products - to be effective as of February 2026*

EP General Chapters

- *Nasal Preparations (Ph.Eur. 12.0/0676)*
- *Preparations for Inhalation (Ph. Eur. 12.0/0671)*

What the Standards Really Mean?



Breaking down key regulatory testing requirements



For generics, passing Q1/Q2 doesn't guarantee equivalent spray characterization performance



Timing of tests in development - why early-phase characterization matters

US FDA Regulatory Requirements

| Drug Product Characterization Studies (During Development) | Specifications for the Drug Product (Release) | Stability |
|--|--|--|
| Priming and Repriming in Various Orientations | Description | Description |
| Effect of Resting Time | Identification | Assay |
| Temperature Cycling | Assay | Impurities and Degradation Products |
| In Vitro Dose Proportionality | Impurities and Degradation Products | Preservatives and Stabilizing Excipients Assay |
| Cleaning Instructions | Preservatives and Stabilizing Excipients Assay | Pump Delivery |
| Device Robustness | Pump Delivery | Spray Content Uniformity (SCU) |
| Effect of Dosing Orientation | Spray Content Uniformity (SCU) | Spray Pattern |
| Profiling of Sprays Near Container Exhaustion | Spray Pattern and Plume Geometry | Droplet Size Distribution |
| Effect of Storage on the Particle Size Distribution | Droplet Size Distribution | Particle Size Distribution (Suspensions) |
| Plume Geometry | Particle Size Distribution (Suspensions) | Particulate Matter |
| Photostability | Particulate Matter | Microbial Limits |
| Preservative Effectiveness and Sterility Maintenance | Microbial Limits | Weight Loss |
| Stability of Primary (Unprotected) Package | Net Content | Leachables |
| | pH | pH |
| | Osmolality | Viscosity |
| | Viscosity | |

EU Regulatory Requirements

https://www.ema.europa.eu/en/documents/scientific-guideline/guideline-pharmaceutical-quality-inhalation-nasal-medicinal-products-revision-1_en.pdf



Table 5.2.1. Pharmaceutical development studies for nasal medicinal products.

| Pharmaceutical development study | Pressurised metered-dose nasal spray | Nasal powders, device-metered | | Nasal liquids | | | |
|--|--------------------------------------|-------------------------------|------------|-------------------|------------------|-------------------|--|
| | | Single-dose | Multi-dose | Single-dose drops | Multidose drops | Single-dose spray | Non-pressurised multidose metered-dose spray |
| (a) Physical characterisation | Yes | Yes | Yes | Yes | Yes | Yes | Yes |
| (b) Minimum fill justification | Yes | No | Yes | Yes | Yes | Yes | Yes |
| (d) Extractables / leachables | Yes | No | No | Yes | Yes | Yes | Yes |
| (f) Particle / droplet size distribution | Yes | Yes | Yes | No | No | Yes | Yes |
| (g) Uniformity of delivered dose | Yes | Yes | Yes | No | No | Yes | Yes |
| (j) Actuator / nasal applicator deposition | Yes | Yes | Yes | No | No | Yes | Yes |
| (l) Shaking requirements | Yes ^a | No | No | Yes ^a | Yes ^a | Yes ^a | Yes ^a |
| (m, n) Initial & re-priming requirements | Yes | No | No | No | No | No | Yes |
| (o) Cleaning requirements | Yes | No | Yes | No | Yes | No | Yes |
| (p) Low temperature performance | Yes | No | No | No | No | No | No |
| (q) Performance after temperature cycling | Yes | No | No | No | No | Yes | Yes |
| (r) Effect of environmental moisture | Yes | Yes | Yes | No | No | No | No |
| (s) Robustness | Yes | Yes | Yes | Yes | Yes | Yes | Yes |
| (t) Delivery device development | Yes | Yes | Yes | Yes | Yes | Yes | Yes |
| (u) Preservative effectiveness / efficacy | No | No | No | No ^b | Yes ^c | No ^b | Yes ^c |
| (x) Spray pattern / plume geometry | Yes | Yes | Yes | No | No | Yes | Yes |

^a For suspensions.

^b Single use formulations should preferably be preservative free, but if a preservative is present, it should be adequately justified.

^c If a preservative is present.

Additional Special Considerations

Unit-dose vs. bi-dose devices:

- Actuation force

Multi-Dose:

- One-time studies - end of life, priming/repriming, in-use study

Component Testing:

- Incoming QC Requirement
 - Challenges with avoiding Release testing even when component QC is robust
- Repeatability Testing

Population Bioequivalence (PBE) versus Average Bioequivalence (ABE)

| ASPECT | PBE | ABE |
|---------------------------|--|---|
| Primary Focus | Compares mean and variability (spread) between Test and Reference. | Compares mean PK parameters (C _{max} , AUC) between Test and Reference. |
| Statistical Basis | Statistical criterion incorporates both mean differences and variance differences. | 90% CI for Test/Reference mean ratio must be within 80-125%. |
| Variability Consideration | Explicitly accounts for differences in variability between products. | Ignores differences in variability between products. |
| When Used | USA | EU |
| Goal | Ensure products have similar average performance and similar distribution of responses. | Ensure products have similar average performance . |

Best Practices in Spray Characterization



Development of spray characterization methods are essential, but complex



“Weight of Evidence” approach: Regulatory expectation for robustness and reproducibility of in vitro methods



Method Validation and Setting Specifications

Spray Actuation Content - Key Considerations



Critical Role of Actuation Parameters

- Stroke length, speed, and consistency directly affect spray output.
- Under- or over-stroking can cause deviations, particularly at end-of-life doses.


Industry Confusion:

- Current **USP <601>** guidance and **Product-Specific Guidances (PSGs)** are not fully aligned.


Method Validation is Essential:

- **Prove your method:** demonstrate reproducibility, accuracy, and robustness.

Acceptable vs. Out of Trend - Spray Pattern



Acceptable vs. Out of Trend - Plume Geometry



Acceptable vs. Out of Trend - DSD



Warning ⚠

If this occurs the measurement data will continue to be analysed and a warning message, denoted by an exclamation mark (!), will be displayed in the Measurement parameters record view.


The reason for the warning will be displayed at the bottom of the measurement parameters window.

Error ⚠


If this occurs the measurement data will not be analysed and red bars are displayed in place of the expected results in the Particle Size History view. Messages are also displayed in the Measurement parameters window, where they are denoted by three exclamation marks (!!!).

The reason for the warning will be displayed at the bottom of the measurement parameters window.


Weight of Evidence




Force to Actuate




0%G Spray Pattern & Deposition




Dmax: 34.91mm Ovality: 1.085



| Method 1 - Vel: 70mm/s Accel: 5000mm/s ² | |
|---|------------|
| LOCATION | DEPOSITION |
| VESTIBULE | ✗ |
| TURBINATES | ✗ |
| NASOPHARYNX | |
| RUN OUT NOSTRIL | ✗ |




Dmax: 39.65mm Ovality: 1.058




| Method 2 - Vel: 130mm/s Accel: 8000mm/s ² | |
|--|------------|
| LOCATION | DEPOSITION |
| VESTIBULE | ✗ |
| TURBINATES | ✗ |
| NASOPHARYNX | ✗ |
| RUN OUT NOSTRIL | |


60%G Spray Pattern & Deposition




Dmax: 19.05mm Ovality: 1.108



| Method 1 - Vel: 70mm/s Accel: 5000mm/s ² | |
|---|------------|
| LOCATION | DEPOSITION |
| VESTIBULE | ✗ |
| TURBINATES | ✗ |
| NASOPHARNYX | |
| RUN OUT NOSTRIL | |



Dmax: 26.56mm Ovality: 1.438



| Method 2 - Vel: 130mm/s Accel: 8000mm/s ² | |
|--|------------|
| LOCATION | DEPOSITION |
| VESTIBULE | ✗ |
| TURBINATES | ✗ |
| NASOPHARNYX | |
| RUN OUT NOSTRIL | ✗ |

Validation of Spray Characterization Methods



Robustness - Method performs reliably under small, deliberate variations in conditions.



Qualification - Demonstrates the method is suitable and fit for its intended purpose.



Repeatability - Produces consistent results when tested under the same conditions.



Intermediate Precision - Confirms reproducibility across different days, analysts, or instruments.

From Data to Decisions

Setting Sample Sizes and Specifications for Nasal Spray Products

Understand
Variability

Define Sample Size
& Tiering

Select Statistical
Method

Justify via Risk &
Data


Key Considerations

- Account for inherent product variability
- Plan for occasional outliers (“rogue devices”)
- No single pass/fail rule: regulators expect a science- and risk-based justification
- Tier 1 vs. Tier 2 testing: balance number of devices and replicates per device (especially for multi-dose)
- Critical Quality Attributes (CQAs): Specs should be tied to attributes that impact clinical performance

Options for Specifications

- Tolerance Intervals (to capture a % of future data)
- Percentile-based limits (e.g., 95th % bounds)
- Clinical relevance anchor (link to PK/PD or deposition performance)
- Historical batch data (development + stability justification)
- Risk-based acceptance criteria (tighter for CQAs, wider for non-critical attributes)

Ensuring Success from the Development Lifecycle



Industry Pain Points & Unmet Needs




Gaps between regulatory guidances and real-world feasibility

Guidance's staying in *Draft* for decades - "Living Document" reality

Not all compendial aspects are reflected in FDA PSGs

Nasal Spray Testing

Next-Generation Approaches



Patient-Centered Testing

Leverage human-realistic actuation profiles that reflect true patient use
Expand physiologically relevant nasal cast models for predictive performance

Integrated Strategies

- Combine in vitro, in silico, and in vivo methods for holistic understanding
- Anticipate and align with evolving regulatory standards

Collaborative Innovation

Strengthen partnerships between industry, regulators, and suppliers
Share data/methods to harmonize best practices and accelerate progress

Key Takeaways



START SPRAY
CHARACTERIZATION
EARLY - DON'T WAIT
UNTIL VALIDATION
STAGE



KEEP TESTING
APPROACHES
ADAPTABLE TO
EVOLVING GUIDANCE



VALIDATE CRITICAL
METHODS LIKE SP/PG
FOR REGULATORY
CONFIDENCE



ENGAGE WITH
FORUMS LIKE IPAC-RS
TO HELP SHAPE
STANDARDS



NASAL INNOVATION FORUM

*Advancing Science & Shaping
the Future of Nasal Drug Delivery*



Thank you for attending!

Questions?