

Aptar Pharma

Delivering for Patients: Navigating Nasal Routes for Optimal Impact

Reenal Gandhi

September 2025

Confidential &
Proprietary



From formulation
to patient.



Aptar
pharma

Intranasal Market Development Snapshots

2



1998

First Launch
Imigran



Zomig, Instanyl,
Narcan (2016)

2018



2020

Valtoco, Tosymra, Spravato,
Nayzilam, Baqsimi



Opvee, Kloxxado,
Naloxone, Zavzpret

2023

Intranasal
Systemic Market
Expanding



Relatively IN Products have a faster growth rate. Revenues CAGR for nasal was 29% compared to 2% CAGR for oral & 4% of inj
Total Nasal Rev \$2B
Total Injectable Rev \$41B

It's Not Oral.
It's Not IV.
It's Enbumyst.



From formulation to patient.

Confidential & Proprietary

Enbumyst™ (bumetanide nasal spray) 0.5 mg



Intranasal has rapid onset of action



Administration time reduced

- No mixing or preparation by pharmacist or nurse
- Dosing is less than a minute



Safe administration by HCP

- No needles and associated waste
- Less chance of dosing errors as nasal spray is prepackaged compared to vials




Intranasal pathways for systemic and direct to CNS can improve effectivity and off-target delivery




Alternative treatment locations

- Satellite clinic instead of hospital infusion center
- Hospital to home




Intranasal (IN)
FDA approved
dose vs
intramuscular
naloxone (IM)





The Intranasal Combination Product Trilogy

5



Confidential & Proprietary

Device Selection Process

Which device?


- Dosage form (liquid, powder)
- Concentration/Dose volume
- Rheology



Droplet size by
Malvern Spraytec

Can I spray it?

- Shear forces
- Droplet size
- Plume characteristics




Plume geometry by
SprayVIEW (Proveris Scientific)

Compatibility
Device
Formulation

- Integrity
- Aggregation
- Quantification

HPLC-SEC



Deposition


- Where is the target in the nasal cavity?



PADA
(lung powder - mice)




Fine Mist Sprayer
(lung/nose liquid – mice)



UDS Powder Pediatric
(nose powder – medium/large animals)



VP7-232NE Pediatric
(nose liquid – medium/large animals - multidose)




UDSL Liquid Pediatric
(nose liquid – medium/large animals)

UDSI Marketed Standard for Intranasal Delivery

8

- 30 years from first product launch on the market
- Ready to use, **primeless** device
- Delivery of metered dose
- Dosage of 100 μl with standard filling volume of water
- Reproducible spray performance by **user-independent** energy accumulation

- Address high **reliability** target
- One-hand actuation, can be used in any position
- Aptar **Plunger** compatible with UDSI Design
- Small quantity of compatible Vials could be supplied
- Customized design possible



Confidential & Proprietary

Human factors studies can help understand target patient groups can handle and actuate devices

Help understand how well patients can use devices and satisfy regulatory expectations with regard to HF studies.

Overall results of the manual actuation data captured from the volunteers (n = 49).

	Study Outcome				
Actuation Force (N)	10	26	43	47	50
% who could actuate	100	100	96	88	71
% who found it difficult to actuate	0	0	39	33	20
% men who found it difficult to actuate	0	0	44	32	16
% women who found it difficult to actuate	0	0	33	33	25

Using prototype devices we studied manual actuation of drug delivery system to understand key parameters:

- Volunteers documented difficulties when force was greater than 40N.
- No differences in adults between men and women
- Type of approach can gather useful information on actuation force to understand actuation human factors


Surface
tension


Formulation
viscosity

Device
orifice

Device
trigger force

- Measured trigger force was not impacted by formulation properties or the discharge of liquid through different nozzle orifice diameters.
- Only the mechanical trigger force had a major impact on the maximum force during spray.
- Trigger force is the critical feature to control to ensure users actuate the nasal drug product as intended





Aptar Nasal Drug Delivery Devices



Liquid Platform

Unidose

Systemic delivery

UDS Liquid

BDS Liquid

Targeted delivery

Cerespray

Powder Platform

UDS Powder

Aptar CSP protective packaging

Multidose

Systemic delivery

CPS Family Preservative Free

VP7 + CB18

Targeted delivery

Phase I ready


VP7 + N2B Actuator

Connected Platform


An innovative approach to nasal drug delivery, leveraging digital technologies to enhance efficacy and patient experience, starting in early development phases

Aptar

- Vaccines are usually developed as **parental injections**
 - Adaptation of existing products for nasal delivery
 - Use existing manufacturing & supply chain
- Typically contain **biologics**
 - More sensitive formulations
 - Shelf life / storage issues
- Vaccination programmes require **cost-effective** solutions
 - Reusable multi-dose drug products (vials)
 - Preservative free systems
 - Easy to use for HCP outside of clinical setting
- Targeting **NALT region** at back of nose
 - Immunity is generated / immune response triggered



Nasal Cavity






AeroNose™ Nasal Cast to Quantify Deposition

Goal:

- Determination of % deposition within anterior area

Description of the Nasal Cast :


- 3D representation of the nasal cavity
- Visualization of the spray deposition on each area
- Quantification of the dose on each nasal area



Block 1 : Nose
Block 2 : frontal sinus ; nasal valve
Block 3 : maxillary sinus ; frontal sinus ; floor of nasal cavity ; turbinates ; ethmoids
Block 4 : maxillary sinus ; sphenoids ; floor of nasal cast ; ethmoids ; turbinates
Block 5 : sphenoids ; rhino-pharynx
Filter : lung fraction (particles < 10µm)



- Commercialized device, which has shown consistent deliver to olfactory Ready to use/Primeless
- Active device (patient does not inhale) ensures accurate, consistent dose
- Powder properties influence the spray performance and regional deposition
- Consistent deposition despite changes to orientation



30°	45°	60°	
28%	30%	35%	0°
35%	48%	46%	5°
	26%		10°

% deposition in olfactory region at difference device orientations (10mg)

UDS Powder Disposable



Max. fill volume: 130 mm³
Fill weight: 10 - 80 mg
Emitted dose weight: >90%

Consider the Patient

Device Selection Studies – An input Formative HF studies, IFU development, Risk Management, etc.




Multidisciplinary approach linking user experience research, human factors regulatory studies and onboarding solutions, accelerating and derisking combination product development for our partners



Noble Parkinson & Alzheimer's User Study

- Noble, an Aptar company specializes in offerings to improve patient lives and support the commercial success of brands in the market
 - *Training devices* – Building patient confidence and preparedness through hands on practice
 - *Services* – Developing user centered solutions through user experience, human factors and design
 - *Solutions* – Driving engagement and access to Noble's offerings through channel partners



17




CONFIDENTIAL | © NOI

Study Structure



- Introduction
- Patient/Caregiver questions
- Hand size measurement
- Lateral pinch strength measurement
- Product familiarization
- Product simulation for measurements – Patient/Caregiver
- UX feedback for both the products

Proprietary




- Dosimetry study in NHPs with [68Ga]NOTA-insulin utilizing an Aptar CPS delivery system to quantify insulin distribution
- Intranasal [68Ga]NOTA-insulin was administered to anesthetized healthy adult vervet monkeys
- Brain regional activity and whole-body dosimetry following PET/CT scans were measured in vervet monkeys
- Aptar CPS successfully delivers [68Ga]NOTA-insulin to the brain
 - Demonstrates regional transport via the olfactory and trigeminal nerves into the brain
- Dosimetry data is favorable for planned future studies of intranasal [68Ga]NOTA-insulin administration in humans

Solingapuram Sai, Craft et al., First biodistribution study of [68Ga]NOTA-insulin following intranasal administration in adult vervet monkeys, resented at AAIC 24




NHP positioning during (1) and after administration (2). Representative a) axial, b) sagittal, and c) coronal brain images after intranasal administration of [68Ga]NOTA-insulin. Red represents highest uptake values.



Enable the Patient to Start Sooner and Stay Longer on Therapy

19

Our state-of-the-art products, development services and patient experience expertise enable pharma to accelerate their go-to-market timelines while promoting improved patient outcomes



COMPREHENSIVE SET OF PRODUCTS, SERVICES AND DATA

ENHANCED PATIENT EXPERIENCE

Journey mapping, educational material, human factors studies, patient value stories...

SERVICES

Digitally-enabled solutions, RWD/RWE generation...


Training and onboarding devices

PRODUCTS

DATA



- Personalized treatment
- Training to technique
- Medication adherence and monitoring
- Behavioral/lifestyle coaching
- Seamless access and support



- Nasal mucosa offers an alternative pathway to the brain that bypasses the BBB -> Targeted delivery to site of action
- Potential to open new therapeutic areas and new molecules to nasal administration
- Highly attractive therapeutic areas fall within the CNS (Central Nervous System) drug delivery arena: Neurodegenerative diseases or Psychiatric conditions ie. Alzheimer's, Parkinson's, depression
- Biologics and complex molecules are showing early promising results with antivirals, mAbs, and peptides. -> Wake Forest Insulin publication



Evolution in systemic and CNS drug delivery—
opportunity to treat unmet needs

Evidence of direct nose to brain transport

Supporting evidence for direct nose to transport in man

Need to further understand transport mechanisms

Opportunities for self administration of vaccines

Incorporate the patient into drug development journey